

Spectroscopic and Electrochemical Investigations of N-(Phosphonomethyl)glycine (glyphosate) and (Aminomethyl)phosphonic Acid (AMPA)

A. Habekost*

Department of Chemistry, University of Education, Reuteallee 46, D-71634 Ludwigsburg, Germany

*Corresponding author: a.habekost@t-online.de

Abstract N-(Phosphonomethyl)glycine (glyphosate), known by the trade name of Roundup[®], is a broad-spectrum systemic herbicide used to kill several types of grass weed. It was first synthesized in 1970 by J. E. Franz, a chemist at the agrochemical corporation Monsanto. Glyphosate's mode of action is to inhibit a plant enzyme involved in the synthesis of some aromatic amino acids. There is some controversy at present about the use of Roundup[®] because its hazard potential is not clear. In this article, we present some reliable and easily performed spectroscopic and electrochemical measurements to identify glyphosate isolated as well as in some commercial products. The analogous experiments apply to (Aminomethyl)phosphonic acid (AMPA), the hydrolysis product of glyphosate.

Keywords: Four year undergraduate/beginners of PhD students, analytical, electrochemistry, IR, MS, UVVISpectroscopy, hands-on learning/manipulatives

Cite This Article: A. Habekost, "Spectroscopic and Electrochemical Investigations of N-(Phosphonomethyl)glycine (glyphosate) and (Aminomethyl)phosphonic Acid (AMPA)." *World Journal of Chemical Education*, vol. 3, no. 6 (2015): 134-140. doi: 10.12691/wjce-3-6-2.

1. Introduction

Glyphosate is one of the most commonly employed herbicides, used on several types of plants, but its hazard potential is still not clear. It is far beyond the scope of this article to report the various arguments for and against the use of glyphosate.

Current methods for analytical measurement of glyphosate and its hydrolysis product AMPA are rather lengthy, expensive, and elaborate, e.g. HPLC, electrophoresis [1] and GC-MSD [2,3]. The latter method requires lowering of the polarity and therefore enhancement of the volatility, which may be done by derivatization of glyphosate and AMPA [2,3], as will be described later.

Kodama et al. [1] used capillary electrophoresis to investigate the contamination of various tea-based beverages with glyphosate. They created a complex of glyphosate with Cu(II) during the electrophoretic run and measured the retention time as a function of pH.




Chiu et al. [4] used capillary electrophoresis as well, but with electrochemiluminescence detection, for the analysis of both glyphosate and AMPA. The detector was an indium-tin oxide (ITO) working electrode in an alkaline phosphate buffer containing Ru(bpy)₃²⁺. In the presence of glyphosate, the anodic wave of Ru(bpy)₃²⁺ in cyclic voltammetry (CV) increased dramatically while the cathodic one decreased. The luminescence signal increased simultaneously.

Sierra et al. [5] studied the electrooxidation of glyphosate on nickel and copper surfaces. The anodic



current peaks in CV increased with the glyphosate concentration. The authors suggested that the detection limit for copper is much lower than for nickel.

Sheals et al. [6,7] successively protonated glyphosate and compared the IR spectra and the extended X-ray absorption (EXAFS) with ab initio calculations. They postulated that the copper ion lies in the center of a (Jahn-Teller distorted) octahedron with all three donor groups of glyphosate. The amine, carboxylate, and phosphonate ligand form two chelate rings with five members each in the equatorial plane.

Three pK_a values for glyphosate are postulated [8]:



The analog values for AMPA are [8]:



Piccolo et al. [9] described the IR spectra of glyphosate at different pH levels, and suggested that the IR bands vary with the dissociation of the glyphosate acid group.

The authors specifically monitored the P-O stretching of the P-OH and P=O groups of the phosphono radical with IR spectroscopy. The given wave numbers have been helpful in assigning the different bands in our own experiments.

Daniele et al. [10] published thermodynamic (stability constants, heat of formation, entropies of formation, and free enthalpies of formation) and spectrophotometric data (λ_{\max} values) of copper (II) complexes of glyphosate in an aqueous solution at different pH levels and temperatures. The data can be explained by assuming the following pH-dependent species: CuLH , CuL^- , CuLH_1^{-2-} , CuL_2^{4-} , and Cu_2L^+ , where L indicates the glyphosate ligand. The authors saw a slight downward shift in λ_{\max} of about 30 nm with increasing pH.

Unfortunately, glyphosate is not electroactive at accessible potentials. Therefore, Pintado et al. [11] electrodeposited copper onto a carbon electrode to determine the concentration of glyphosate in drinking water. They saw an enhancement of the current in CV on adding glyphosate solution, due to the formation of a copper-glyphosate complex. The aim of their paper was "to present an electrochemical method for the quantification of glyphosate in a simple, rapid, and inexpensive way." By adding glyphosate to the copper glassy-carbon electrode, the anodic current peak increased linearly with the glyphosate concentration. The copper on the carbon-paste electrode showed a different behavior: Two oxidation peaks were found, which increased with increasing glyphosate concentration due to the formation of Cu^+ - and Cu^{2+} -glyphosate complexes. The cathodic peaks also increased, and simultaneously the peaks shifted to more negative potentials on increasing the glyphosate concentration.

Dos Santos et al. [12] used square wave voltammetry (SWV) of a hanging mercury-drop electrode to determine the concentration of glyphosate at different levels of pH. The glyphosate formed complexes with copper ions. The current peak in CV was measured as a function of the copper (II) concentration. In addition, they measured the current as a function of the deposition time of the Cu(II)-glyphosate complex adsorbed onto the electrode.

Börjesson and Torstensson [2], following Deyrup et al. [3], developed new methods to determine the concentration of glyphosate and AMPA with GCMS, such as derivatizing both molecules with a mixture of trifluoroacetic anhydride and trifluoroethanol (volume 2:1). Afterwards, the solution was kept at 100°C for an hour. After cooling to room temperature, the sample was evaporated, redissolved in 1 mL ethylacetate, and then analyzed with GC-MSD. This method is described in details later in this paper.

With such experiments, Börjesson and Torstensson [2] could quantify glyphosate and AMPA near a railway line in Sweden more than a year after exposing the area. Furthermore, the authors estimated the limit of detection of glyphosate in the soil to be approximately 0.003 mg/g.

At first glance, it does not seem to be easy to determine the concentration of glyphosate and its main metabolite AMPA in a simple way, in particular for inexperienced students. But the derivatization method, the formation of copper complexes, and the detection of glyphosate and AMPA with GC-MSD, CV, and UVVIS spectroscopy are practicable ways to improve the students' analytical knowledge. In our experience, advanced students can carry out all of these experiments in about 40 hours.

First, we describe the pH-dependence of glyphosate and AMPA by titration and determine the pK values. Then, we describe the production and the investigation of the

colored Cu(II)-glyphosate complexes by UVVIS spectroscopy at various levels of pH. These complexes can also be simply identified by IR spectroscopy. We then experimentally determine glyphosate and AMPA with GC-MSD after derivatization with TFAA and TFE. Finally, following Sierra [6,7] and Börjesson [2], we present the results of electrochemical quantification of glyphosate and AMPA as well as for glyphosate in Roundup®.

2. Pedagogical Objectives

Chemistry students need thorough training in analytical theory and practice to improve their chemical knowledge. Therefore, the aim of this article is a versatile experimental investigation of two important environmental chemicals. This includes modern methods such as spectroscopic and electrochemical knowledge in acquisition and interpretation of mass spectra, IR spectra, and cyclic voltammograms. UVVIS spectroscopy and titration are more classical but are still very useful methods. These methods, easy to use and quick to carry out, are the first step in collecting analytical information.

With respect to environmental analysis, the quality of the applied analytical method for the substances under investigation is determined by the recovery rate of the substances in their natural environment. Therefore, one has to first extract the substance, which is normally adsorbed in soil, by the Soxhlet extraction method before the analytical methods can be applied.

We think that the pedagogical benefit of the described procedures is that investigation of the same substances with different analytical methods is enabled. Furthermore, a comparison between the pure substance and the complex mixture of the commercial products may be instructive.

3. Experimental Procedures

Chemicals and instruments:

Glyphosate (Sigma Aldrich, 54521)
 AMPA (Sigma Aldrich, 324817)
 Potentiostat (μ -Stat 400, DropSens)
 Screen-printed electrodes (DS 550)
 KOH solution (10%)
 Sodium carbonate (saturated)
 Double-distilled water
 Copper sheets or CuO powder (Hedinger, Germany)
 Bottles with snap-on caps
 Photometer (Perkin Elmer Lambda XLS+)
 FTIR spectrometer (Bruker, Vector 22)
 KBr pellets

3.1. Preparation of Cu-glyphosate and Cu-AMPA

The aim of our experiments in general is to simplify the experimental approach for the detection of glyphosate and AMPA, rather than obtaining preferably low detection limits.

Therefore, we took conventional copper sheets (0.5 g \approx 0.0077 mol) or CuO powder (0.62 g \approx 0.0077 mol), and put them into solutions with different pH levels (5g each

of water and sodium carbonate). We used the following chemical approaches:

1. Cu in water
2. Cu in sodium carbonate
3. CuO in water
4. CuO in sodium carbonate

Solutions both with and without glyphosate (0.02 g \approx 0.00012 mol) or AMPA (0.015 g \approx 0.00014 mol) were used. Quantity of copper was significantly higher compared to glyphosate or AMPA.

Figure 1 shows pictures of Cu in water with glyphosate and AMPA respectively, Cu in sodium carbonate with glyphosate and AMPA respectively, and Cu only in sodium carbonate (right) after 12 hours.

It is obvious that only if glyphosate or AMPA are present, the color changes to blue. The color is deeper in sodium carbonate than in water and with glyphosate than with AMPA.



Figure 1. left: Copper sheets in water with glyphosate (left) and with AMPA (right); Right: Copper sheets in sodium carbonate with glyphosate (left), with AMPA (right), and without glyphosate or AMPA (middle)




Figure 2. Spectra of Cu-glyphosate and Cu-AMPA




Figure 3. CuO in sodium carbonate, Cu in sodium carbonate, Cu in KOH, and Cu in H₂O after 12 hours, respectively

In Figure 2, the spectra of glyphosate and AMPA in water are shown. The spectra match those of Daniele and Coutinho [11,13]. This means that λ_{\max} (Cu-AMPA) \approx 680 nm, and λ_{\max} (Cu-glyphosate) \approx 720 nm.

Figure 3 demonstrates the above-mentioned shift of λ_{\max} on increasing the pH by about 30–40 nm (comparing H₂O and KOH as solvents).

Furthermore, the formation of the Cu-glyphosate complex is more pronounced in sodium carbonate than in water.

3.2. Titration

Chemicals and instruments:

Sodium hydroxide (0.1 mol/l)

pH electrode

Dosimeter (Metrohm, 554)

Data acquisition system (Sensor Cassy, LeyboldDidactic, Germany)

Stirrer

Burette

In the following experiment, we attempt to determine the pK_a values of glyphosate and AMPA (Figure 4). First, we acidify the aqueous solution with 2M hydrochloric acid and titrate with 0.1 M sodium hydroxide solution.

The experimental pK_a values comply with the values mentioned above. The corresponding pK_a values of AMPA are smaller than those of glyphosate.




Figure 4. Titration of glyphosate and AMPA

Therefore, one can titrate some commercial products to quickly test whether they contain glyphosate or not. In Figure 5, one can clearly see that Weed killer[®] and Roundup[®], two tested products, contain glyphosate.




Figure 5. Titration of Weed killer[®] and Roundup[®] compared to pure glyphosate

3.3. FTIR Spectra

Figure 6 shows the FTIR spectrum of glyphosate and AMPA in a KBr pellet (300 mg KBr and 1 mg substance; pressure: 1.5 t).




Figure 6. FTIR of glyphosate (top) and AMPA (bottom) in a KBr pellet.

In Table 1a and Table 1b, we report the assignments of the various bands of glyphosate and AMPA.

Table 1a. IR bands assignments (cm^{-1}) for glyphosate between 2000 and 900 cm^{-1} [6,7]

| | | | |
|------------------------------|------|--------------------------|------|
| C=O of free COOH | 1732 | CH ₂ groups | 1243 |
| C=O of H-bonded COOH | 1720 | P-OH | 1223 |
| C-O asymmetric vibration | 1568 | | |
| NH ₂ deformations | 1557 | CO, OH groups | 1203 |
| NH ₂ deformations | 1483 | P-OH | 1170 |
| CO and OH of free COOH | 1434 | P-O ⁻ | 1090 |
| CO and OH of H-bonded COOH | 1424 | CCNC skeletal vibrations | 1081 |
| CH ₂ deformations | 1336 | CCNC skeletal vibrations | 1031 |
| P-O of PO ₃ H | 1268 | P-OH | 1000 |
| | | CCNC | 916 |

Table 1b. Some IR bands assignments (cm^{-1}) for AMPA between 2000 and 400 cm^{-1} [12]

| | |
|-------------------------------|------------------|
| NH ₂ deformations | 1650, 1620, 1530 |
| CH deformations | 1443 |
| P=O | 1164 |
| PO ₃ | 1031 |
| P-C | 727 |
| HO-P=O and O=P=O deformations | 464 / 452 |

The carboxylate band of glyphosate is detected roughly between 1300 and 1700 cm^{-1} , whereas the bands between 950 and 1200 cm^{-1} and around 1600 cm^{-1} originate from the phosphonate and the amine group respectively. Complexation of glyphosate can thus be studied by the analysis of these bands. (Note that the glyphosate bands at 1732 cm^{-1} and 1424 cm^{-1} represent the undissociated COOH group).




Figure 7. FTIR of glyphosate (black), Cu-glyphosate after 10 minutes (red), and Cu-glyphosate (green line) after 12 hours

Figure 7 shows the increasing complexation of glyphosate with Cu(II). After 10 minutes (red curve), the C=O vibration of the COOH group splits to 1602 cm^{-1} and the band at 1424 cm^{-1} shifts to about 1400 cm^{-1} , according to Sheals [6,7]. At the end of the complexation, the shift to lower frequencies is completed. Furthermore, the phosphonate and the amine groups are bonded to Cu(II) as well, causing a shift in the P-O bands (around 1100 cm^{-1}), while the CH deformations (1336 cm^{-1}) are not affected.




Figure 8. FTIR spectrum of AMPA (black) and Cu-AMPA (red)

The comparison between AMPA (black line in Figure 8) and Cu-AMPA (red line in Figure 8) shows a shift of the NH₂ deformation band to lower frequency numbers. The P-O frequency numbers (e.g. the PO₃ band) vary in the lower frequencies as well, as expected.

3.4. GC-MSD Measurements

Chemicals and materials:

- Trifluoroacetic anhydride (SigmaAldrich, No. 106232)
- Trifluoroethanol (SigmaAldrich, No. T 63002)

GC-MSD with autosampler (GC: Hewlett Packard 5890, MSD: Hewlett Packard 5972, autosampler: Hewlett Packard 6890)

Column: RTX-35

Carrier gas: He 5.0

50-mL round-bottom flask

Reflux condenser

Capillary air bleed for solvent evaporation.

Experimental procedure:

10 μg glyphosate was mixed with 200 μg trifluoroacetic anhydride (TFAA) and 100 μg trifluoroethanol (TFE) in a 50-mL round-bottom flask. The mixture was heated to 90°C with a controlled heater and refluxed for about an hour.

Afterwards, the mixture was flushed with clean air for two minutes to evaporate the solvents. The residue was then mixed with 10 mL ethyl acetate and directly analyzed with GC-MSD.

The temperature profile of the GC was 100°C in the first minute and a final temperature of 200°C. Temperature increased at a rate of 5°C/minute at an injection volume of 1 μL . The MS detector started after the solvent peak at 2.5 minutes.

The GC-MS detection was performed according to Deyrup et al. [3] and Börjesson et al. [2]. As described above, both research groups derivatized the trifunctional

glyphosate and the bifunctional AMPA by simultaneous esterification and acylation of the carboxylate, the phosphonic, and the amino group.

In the GC of glyphosate (Figure 9), one can see a remarkable peak at 8 minutes, and two other peaks at 5.8 minutes and 7.2 minutes. The first peak corresponds to a contamination of the injection syringe, while the second peak was formed either due to the production of the glyphosate derivative or the thermal destruction during the GC analysis. The EI-mass spectrum of the peak at 7.2 minutes (molecular peak m/z 464) may be due to the loss of a P-O group

The corresponding EI-mass spectrum of glyphosate is in accordance to those described in Börjesson and Deyrup [2,3]: The molecular ions of the derivatives of glyphosate (m/z 511) show only low intensities, whereas the peak m/z 411 (loss of TFE) has a strong intensity. Other fragments are: m/z 492 (loss of HF), m/z 411 which is indicative for $\text{C}(\text{O})\text{CH}_2\text{N}[\text{C}(\text{O})\text{CF}_3]\text{CH}_2\text{P}(\text{O})(\text{OCH}_2\text{CF}_3)_2$, and m/z 113 which is an unidentifiable fragment.

Similarly, the molecular peak of the AMPA derivative (m/z 371) is not strongly developed. The fragment peaks (m/e 302: loss of CF_3 ; m/e 246: $(\text{CF}_3\text{CH}_2\text{O})_2\text{P-OH}$; and m/e 126: $\text{CF}_3\text{C}(\text{O})\text{NHCH}_3$ fragment) are much more intense.




Figure 9. I: GC spectra of the glyphosate derivative C, A: AMPA (contamination of the injection syringe), and B: glyphosate fragment m/z 464; II: GC and EI-mass spectrum of the AMPA derivative; III: GC-MS spectrum of both AMPA and glyphosate derivatives simultaneously




Figure 10. GC-MS spectrum of Roundup®—glyphosate peak C at about eight minutes.

Figure 10 indicates the difference between pure glyphosate (Figure 9, I) and the commercial, glyphosate-

containing product called Roundup®. We followed the same processing procedure as mentioned above. In case of Roundup®, the GC shows several peaks with retention time above eight minutes (glyphosate), but the glyphosate peak can clearly be identified by its retention time of eight minutes.

3.5. Electrochemical Measurements

Chemicals and materials:

Phosphate buffer (pH 6.9, 0.1 mol)

$\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ solution (10 mmol)

Glyphosate (1 mg)

Potentiostat (μStat 400)

Screen-printed electrodes (DS 550, working electrode:

Pt; counter electrode: Pt; reference electrode: Ag)

Pipettes

Double-distilled water for purification of the electrodes.

The experimental procedure is the same as given in Pintado [11]. First, we dropped 1 mL of a copper sulfate solution (12.5 mg copper sulfate in 5 mL of a 0.1 mol/l phosphate buffer, pH = 6.9) onto the screen-printed electrodes (DS 550: working electrode: Pt; counter electrode: Pt; reference electrode: Ag). Then, we electrodeposited copper onto the platinum working electrode for 60 seconds at -0.8 V. This was carried out by constant stirring. Afterwards, we took a cyclic voltammogram nine times in a row (scan range: -0.8 V \rightarrow 0.6 V \rightarrow -0.8 V; scan rate: 0.1 V/s; no stirring) as shown by the red dotted lines in Figure 11. One can see that the anodic current peak (oxidation of copper to Cu(II)) decreases due to the gradual reduction of the copper layer.

We then mixed 1 mL of the copper sulfate solution with 1 mg glyphosate and executed the same procedure. It is obvious that the current peaks drift to lower potentials (of about 0.12–0.2 V) according to Pintado [11], and increase in intensity (compared to the corresponding CVs without glyphosate). We did not attempt to quantify the detection limit of glyphosate, because we only wanted to test the possibility of easy electrochemical detection of glyphosate.




Figure 11. Cyclic voltammograms of copper sulfate after deposition of copper (red, dotted curves) and of copper sulfate / glyphosate mixture (same conditions). Only the anodic curves are shown. The arrows indicate increase of the CV scan.

3.6. Glyphosate from Soil

Chemicals and materials:

Glyphosate
Double-distilled water
Garden soil
Soxhlet extractor
Round-bottom flask
Heater
Reflux cooler
Rotary evaporator

Experimental procedure:

We dissolved 20 mg glyphosate in 5 mL double-distilled water. Of this, 2.5 mL was removed and dried in a drying oven at 70°C . The residue was derivatized as described above and solved in 1 g ethyl acetate. GCMS was recorded as the “reference.” At the same time, the other 2.5 mL solution was mixed with 15 g soil and extracted with 85 mL water in a Soxhlet extractor for two hours. Then, the extract was concentrated with a rotary evaporator, dried in a drying oven at 70°C , and analyzed with GCMS after derivatization, being the “probe.”

Figure 12 shows the extraction equipment:



Figure 12. Extraction equipment: heater, flask, Soxhlet with extraction thimble containing glyphosate and soil, and the reflux cooler.

The quotient of the integrals of the glyphosate GC peaks with and without soil is: $I_{\text{probe}}/I_{\text{reference}} \approx 91\%$ recovery rate – quite an acceptable result.

4. Conclusion

In this paper we describe several spectroscopic and electrochemical methods to identify glyphosate and AMPA as pure substances, in commercial products and in soil. We restrict ourselves to methods that

- show clear results, and
- are easily to perform.

All experiments can be done in about a week (40 hours) in an analytical lab course. With this plethora of experiments, students not only improve their experimental skill, but also learn fundamental theoretical aspects about UVVIS, IR, GCMS, and cyclic voltammetry.

Therefore, we think that the versatile investigation of glyphosate and AMPA is a significant contribution to professional training of students in analytical chemistry.

References

- [1] Kodama, S., Ito, Y., Taga, A., Nomura, Y., Yamamoto, A., Chinaka, S., Suzuki, K., Yamashita, T., Kemmei, T., and Hayakawa, K, A fast and simple analysis of glyphosate in tea beverage by capillary electrophoresis with on-line copper (II)-glyphosate complex formation, *J. Health Sci.*, 54, 602-606 (2008).
- [2] Börjesson, E. and Torstenson, L, New methods for determination of glyphosate and (aminomethyl) phosphonic acid in water and soil, *J. Chromatogr.*, 886, 207-216 (2000).
- [3] Deyrup, C.L., Chang S.M., Weitraub, R.A., Moye, H.A., Simultaneous esterification and acylation of pesticides for analysis by gas chromatography. 1. Derivatization of glyphosate and (aminomethyl)phosphonic acid with fluorinated alcohols-perfluorinated anhydrides, *J. Agric. Food Chem.* 33, 944-947 (1985).
- [4] Chiu, H.Y., Lin, Z.Y., Tu, H.L., and Whang, C.W, Analysis of glyphosate and aminomethylphosphonic acid by capillary electrophoresis with electroluminescence detection, *J. Chromatogr A*, 1177, 195-198 (2008).

- [5] Sierra, E., Mendez, M.A., Sarria, V.M., and Cortes, M.T. Electrooxidación de glifosatosobre electrodos de níquel y cobre, *Quim. Nova*, 31, 220-226 (2008).
- [6] Sheals, J., Persoon, P., and Hedman, B., IR and EXAFS spectroscopic studies of glyphosate protonation and copper (II) complexes of glyphosate in an aqueous solution, *Inorg. Chem.*, 40, 4302-4309 (2001).
- [7] Sheals, J., Granström, M., Sjöberg, S., and Persson, P., Coadsorption of Cu(II) and glyphosate at the water-goethite (α -FeOOH) interface: molecular structures from FTIR and EXAFS measurements, *J. Colloid and Interface Science* 262, 38-47 (2003).
- [8] http://www.eurl-pesticides.eu/library/docs/srm/1stws2006_lecture_christensen_glyphosate.pdf. Delminoca, E.L., Bertozzi, J., de Souza, N.E., and Oloveira, C.C., Determination of glyphosate and aminomethylphosphonic acid for assessing the quality tap water using SPE and HPLC, *Acta Scientiarum*, 36, 513-519 (2014).
- [9] Piccolo, A. and Celano, G., Modification of infrared spectra of the herbicide glyphosate induced by pH variation, *J. Environ. Sci. Health*, B 28, 447-457 (1993).
- [10] Daniele, P.G., De Stefano, C., Prenesti, E., and Sammartano, S. Copper (II) complexes of N-(phosphonomethyl)glycine in aqueous solution: a thermodynamic and spectrophotometric study, *Talanta*, 45, 425-431 (1997).
- [11] Pintado, S., Amaro, R.R., Mayen, M., and Mellado, J.M.R. Electrochemical determination of the glyphosate metabolite aminomethylphosphonic acid (AMPA) in drinking water with an electrodeposited copper electrode, *Int. J. Electrochem. Sci.* 7, 305-312 (2012).
- [12] Dos Santos, S.C., Galli, A., Felsner, M.L., and Justi, K.C., Desenvolvimento de metodologia eletroanalítica para a determinação do pesticida glifosato em amostras ambientais, *Revista virtual de Química*, 6, 866-883 (2014).
- [13] Coutinho, C.F.B., Coutinho, L.F.M., Mazo, L.H., Nixdorf, S.L., Camara, C.A.P., and Lancas, F.M., Direct determination of glyphosate using hydrophilic interaction chromatography with coulometric detection at copper microelectrode, *Anal. Chim. Act.*, 592, 30-35 (2007). <http://pubs.usgs.gov/tm/tm5a10/pdf/tm5a10.pdf>.